

NeedleVISE®


5-Port and 1-Port
sharps securement

The Smart Way
to Manage Sharps



NeedleVISE adds safety when the use of integrated safety needles is not possible. The NeedleVISE secures sharps, eliminates tray penetration, provides for easy needle insertion and accepts a wide range of sizes.

NeedleVISE®



NeedleVISE 5-Port at a glance

- Holds five needles 17ga to 27ga
- Ensures needles can't penetrate tray
- Organizes sharps for fast accountability
- Compact design for easy disposal in sharps container
- Fits into your procedure tray



NeedleVISE 1-Port at a glance

- Universal solution for safe needle securement
- Only device that secures large gauge needles
- Holds sharps 8ga to 27ga
- Large target for safe and easy insertion

Atrion[®]
Medical

Atrion Medical Products, Inc. • 1426 Curt Francis Rd • Arab, AL 35016 USA
+1 (256) 317-2233 • info@atrionmedical.com • www.atrionmedical.com

Protected under U.S. Patent 6,382,417. International Patents Pending.
NeedleVISE® is a registered trademark of Atrion Medical Products, Inc.